INDEX

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SURFACES SELECTION GUIDE

surface	code	binding capacity	binding interaction	Biomolecule properties	suggested uses
Medium Binding	MB	100-200 ng IgG/cm ²	A polystyrene surface less polar than the HB 8 one and with affinity to molecules of a more hydrophobic nature. The MB surface can better work with lipidic molecules	Large molecules (> 20 kD) with large or abundant hydrophobic regions	Assays where the molecules that have to be coated present large hydrophobic regions or assays that need hydrophobic polipeptides as coated molecule
High Binding	HB8	Up to 400-500 ng IgG/cm ²	A polystyrene surface with high affinity to molecules with mixed hydrophilic/hydrophobic domains such as proteins and antibodies. The surface is optimized for antibodies (IgG)	Improves binding of medium to large molecules (> 10 kD) that are positively charged with or without hydrophobic regions	Assays where the adsorbed molecule must work in excess (up to 400-500 ng/cm ²) - e.g. ELISAs (indirect ELISA test for IgG or ELISA capture test for IgM detection) Very suitable to set up competitive tests (molecule < 50 ng/cm ²) -e.g. steroids/hormones - for the ability to orientate the adsorbed molecule
NO Binding	NB	N.A.	None-Inhibits hydrophobic and ionic interactions	Significantly reduces protein binding	Certain procedures require a surface that is non binding. Many proteins,enzymes in particular may become activated or inactivated upon attachment to a surface. A non-binding surface can be ideal for assays requiring this feature
Streptavidin coated	SA 5/200	12p/mol of biotin The binding will vary depending on the nature of the biotin-bound moieties and the steric hindrance	Non covalent interaction towards biotinylated molecules; streptavidin-biotin interaction is nearly irreversible since the binding affinity (K _a) in solid-phase is around 10 ⁸ -10 ¹⁰	any biotinylated molecules	Streptavidin coated microplate is designed for analysis of biotinylated molecules or for indirect coating with biotinylated ligands. This support binds very effectively any biotinylated molecule, with very rapid kinetic and excellent binding capacity, offering a multipurpose tool for various capture and immunodetection systems
Protein A coated	PA	~ 5 pmol IgG/well	specific for Fc region of IgG	Most immunoglobulins of mammalian species	- Binds strongly to IgG from human, rabbit,guinea pig,pig, dog and cat - Binds strongly to mouse IgG2a,IgG2b and moderately to IgG3 - Binds to Fc region of antibodies for optimal orientation
Protein G coated	PG	~ 5.3 pmol IgG/well	specific for Fc region of IgG	Most immunoglobulins of mammalian species	- Binds strongly to IgG from human, rabbit,mouse ,pig,bovine,dog , goat and horse - Binds only to IgG- no cross-reactivity with other antibody classes - Binds to Fc region of antibodies for optimal orientation
Poly-L or D-lysine coated	Lys-L Lys-D	N.A.	This polymer creates a uniform net positive charge on the plastic surface	Cells and nucleic acids	This surface can enhance cell attachment, growth and differentiation. Moreover this surface can tight bind negatively charged nucleic acids like dsDNA
Poly-L-Arginine	AR		high density of: α -amino/ α -carboxyl/ guanidino groups able to react trough electrostatic and stereospecific bonds		interactions with serino proteases interactions with maturation promoting factors

biomat

surface	code	binding capacity	binding interaction	Biomolecule properties	suggested uses
Calmodulin	Cal		hydrophobic sites		interactions with proteins involved in glycogen metabolism interactions with factors involved in neurotransmission mechanism interactions with enzymes involved in the NAD*/NADP* phosphorylation system
Biotin	Bio		streptavidin/biotin non covalent biological strongest interaction		interactions with avidin interactions with streptavidin
Concanavalin A coated	Con A	N.A	specific for C3-C4-C5 hydroxyl groups	carbohydrates	High affinity for terminal α -D-mannosyl and α -D-glucosyl residues . Used to immobilise glycoproteins or carbohydrates containing these groups
Jacalin coated	Jac	N.A.	high affinity for α -D-galactose residues of the biochemical structure of IgA1 and cell membranes	Human IgA 1	Human IgA1 specific binding, sterically oriented Purification of human immunoglobulins (especially IgA1) Separation of immunocomplexes antigen-antibody Separation of IgA1 from contaminants Stimulation of T-cells
Wheat germ coated	WG	N.A.	specific binding of N-acetyl-D-glucosamine residue	Glicoproteins, enzymes and cell membranes	Studies of surfaces of normal and transformed cells Glycoprotein purification including membrane glycoproteins Studies of cell surface changes during development and the cell cycle
Aminated with primary and secondary amines	AM 1 AM2	N.A	Reactive amino groups allow immobilization of molecules through reductive amination of aldehyde or ketone groups without prior activation of the surface. Alternatively, molecules can be immobilized through amide-bond formation with carbodiimide-activated carboxyl groups. Bifunctional crosslinkers may be used to introduce other functional groups.	Proteins, peptides, carbohydrates, enzymes	Immobilisation of molecules which are bound weakly or not at all by physical adsorption, namely small peptides (M.W. 1000-5000) drugs, toxins or hormones. Oriented immobilisation of molecules in order to secure the integrity and accessibility of their specific sites avoiding the risk of inhibition of these sites by casual physical adsorption for such molecules as Fab-SH-antibody fragments, or nucleic acids (single or double strand)
Carboxylated	COOH	N.A	Carboxyl groups can be performed with a carbodiimide followed by coupling of an amino containing molecule, resulting in a stable amide bond between the surface and the molecule. Alternatively bifunctional crosslinkers may be used to introduce other functional groups.	Molecules containing amino groups	The amino group presents in any molecules, such as peptides or proteins, binds to COOH surface through formation of amide bonds between the amino group presents in the molecule and the surface carboxylic group by the action of carbodiimide .

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ANALYSIS OF DATA

The data obtained from the two types of samples of microplates were processed in the following way:

NEGATIVE SAMPLES	POSITIVE SAMPLES
a. min. O.D. observed for each type of strip	a. coefficient of correlation
b. Max. O.D. observed for each type of strip	b. linear regression calculated as $y = a + bx \quad (1)$
c. mean of Standard Deviations	with a confidence level of 95% X values were those obtained with competitor's samples Y values were those obtained with <i>biomat's</i> HB8 samples

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Method of analysis

Coating solution: 2.5 µg/ml of Rabbit IgG in 0.1 M Carbonate Buffer pH 9.6.
Conjugate: Goat Anti-IgG-HRP

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approvato
approvato
approvato

Test

visual examination
general appearance
packaging

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Two types of tests were made for checking the properties of our **No Binding** capacity surface :

a. No binding capacity against IgG (met. 14)

1. coating: dispense 100µl/well of Rabbit Anti-HIgM (concentration of 25 µg/ml in Carbonate Buffer pH 9.6) and incubate 24 h at 4°C.
2. 4 washing with PBS pH 7.2
3. blocking: dispense 150 µl/well of BSA 1% in PBS pH 7.2 and incubate overnight at 4°C.
4. wash 4 times with 0.1M PBS pH 7.2 + 0.05% Tween® 20
5. dispense 100 µl/well of Biotinylated Goat Anti- Rabbit IgG (concentration from 0.062 to 0.0007 µg/ml) and incubate 1 hr at room temperature
6. wash 4 times with 0.1M PBS pH 7.2 + 0.05% Tween® 20
7. dispense 100 µl/well of Streptavidin-POD (concentration 0.16 µg/ml) and incubate 1 hr at room temperature
8. wash 4 times with 0.1M PBS pH 7.2 + 0.05% Tween® 20
9. dispense 100 µl/well of TMB
10. after 30' stop the reaction with H₂SO₄ 1 N (100 µl/well)
11. reading is made at 450 nm after 5'

results

O.D. values	NO binding	Medium binding	High binding HB 8
<i>RT</i>	62	974	1563
<i>45° C</i>	60	n.a.	1476
<i>- 20°C</i>	61	n.a.	1590

The test was performed after 7 days in the above conditions

b. No binding capacity against Albumin (met. 25)

1. coating: dispense 100µl/well of a mixture of 20 ng/ml Biotinylated Albumin diluted in 3 µg/ml of BSA in 0.1 M PBS pH 7.2 and incubate 24 h at 4°C.
2. wash 4 times with 0.1M PBS pH 7.2 + 0.05% Tween☒ 20
3. dispense 100 µl/well of 0.1 µg/ml Streptavidin-POD and incubate 30' at room temperature
4. wash 4 times with 0.1M PBS pH 7.2 + 0.05% Tween☒ 20
5. dispense 100 µl/well of TMB
6. after 10' stop the reaction with H₂SO₄ 1 N (100 µl/well)
7. reading is made at 450 nm

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

B. Sulfosuccinimidyl maleimidomethyl cyclohexane carboxylate (SMCC),

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
biomat

C. N-Hydroxysulfosuccinimide (Sulfo-NHS) 또는 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide Hydrochloride(EDC)와 합쳐진 N-Hydroxysuccinimide (NHS).

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1. peptide
2. Sulfo NHS (or NHS)
3. EDC
4. intermediate active compound resulting from the reaction
5. *biomat* NH₂ surface
6. peptide covalently immobilised on *biomat* NH₂ surface

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biomat

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Solid phase:	Biomat plates	MT02F-AM1 (primary amino groups) MT01F-HB (high binding capacity)
d-Biotin	Sigma	Cat. No. B 4501
1-Ethyl-3-(3 dimethylaminopropyl)- carbodiimide (EDC)	Sigma	Cat. No. E 1769
Sulfo-N-hydroxysuccinimide (sulfo-NHS)	Fluka	Cat. No. 56485
Dimethylsulfoxide (DMSO)	Merck	Cat. No. 2931
Tween® 20	Merck	Cat. No. 822184
Streptavidin	BIO-SPA	Cat. No. S002-60
Streptavidin-peroxidase conjugate	BIO-SPA	Cat. No. SB01-61
BSA	Intergen	Cat. No. 3100
TMB peroxidase substrate	Kirkegard & Perry	Cat. No. 50-76-05

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
PRODUCT CATALOGUE

General information above codes	code
96 WELL PLATE 12 X 8 FRAME	MT
96 WELL PLATE SINGLE WELL HOLDING FRAME	MG
96 WELL SOLID PLATE	MC
TUBES	TS

MATERIAL **TRANSPARENT POLYSTYRENE Unless otherwise indicated**
 Strips and plates in black or white polystyrene are available upon request

pack **case**
 sleeve carton

TISSUE CULTURE PRODUCTS

Code	Description	quantity per pack	quantity per case
MC024FL-TC/S	24 well plate Tissue Culture treated- sterile with lid	1	80
MC024F-TC	24 well plate Tissue Culture treated	20	160
MC0FL-TC/S	96 well plate Tissue Culture treated- sterile with lid	1	160
MC0F-TC	96 well plate Tissue Culture treated	25	200
PD 60-TC	Ø 60 mm Petri dish Tissue Culture treated	10	500
PD 90-TC	Ø 90 mm Petri dish Tissue Culture treated	20	500
PD 120-TC	Ø 120 Petri dish Tissue Culture treated	20	360
PD 150-TC	Ø 150 Petri dish Tissue Culture treated	20	160

CRYSTALLIZATION PLATES

Code	Description	quantity per pack	quantity per case
MC0-CR O	96 well protein crystallization plate Under Oil	10	100
MC0-CR 96-2	96 well protein crystallization plate 192 wells	10	100
MC0-CR 96-3	96 well protein crystallization plate 288 wells	10	100

TRANSPARENT BOTTOM PLATES

Code	Description	quantity per pack	quantity per case
MCB-GB-96	96 well plate Black PS - Glass Bottom	5	40
MCB-GB-384	384 well plate Black PS - Glass Bottom	5	40
MCB-GB-1536	1536 well plate Black PS - Glass Bottom	5	40

HTS PLATES

Code	Description	quantity per pack	quantity per case
MCW0F-96	96 well plate White PS	25	200
MCW0F-384	384 well plate White PS	25	200

IMMUNOASSAY PRODUCTS**MEDIUM BINDING CAPACITY**

code	description	quantity per pack	quantity per case
MT0F2-MB	1x8 flat bottom wells strips on 12 x 8 frame	25	200
MG0F-MB	1x8 flat bottom wells strips on single well holding frame	25	200
MC0F-MB	96 wells plate - flat bottom	25	200

HIGH BINDING CAPACITY HB8

code	description	quantity per pack	quantity per case
MT01F2-HB8	1x8 flat bottom wells strips on 12 x 8 frame	25	200
MG01F-HB8	1x8 flat bottom wells strips on single well holding frame	25	200
MC01F-HB8	96 wells plate - flat bottom	25	200

SPECIAL SURFACES PRODUCTS**NO BINDING CAPACITY**

code	description	quantity per pack	quantity per case
MT03F2-NB	1x8 flat bottom wells strips on 12 x 8 frame	25	200
MG03F-NB	1x8 flat bottom wells strips on single well holding frame	25	200
MC03F-NB	96 wells plate - flat bottom	25	200
TS 127503-NB	12 x 75 mm tubes	100	1000-2000

biomat

PRIMARY AMINO GROUPS

code	description	quantity per pack	quantity per case
MT02F2-AM1	1x8 flat bottom wells strips on 12 x 8 frame	10or25	100 or 200
MG02F-AM1	1x8 flat bottom wells strips on single well holding frame	10or25	100 or 200
MC02F-AM1	96 wells plate - flat bottom	10or25	100 or 200
TS 127502-AM1	12 x 75 mm tubes	100	1000-2000

SECONDARY AMINO GROUPS

MT02F2-AM2	1x8 flat bottom wells strips on 12 x 8 frame	10or25	100 or 200
MG02F-AM2	1x8 flat bottom wells strips on single well holding frame	10or25	100 or 200
MC02F-AM2	96 wells plate - flat bottom	10or25	100 or 200
TS 127502-AM2	12 x 75 mm tubes	100	1000-2000

CARBOXYLATED

code	description	quantity per pack	quantity per case
MT04F2-COOH	1x8 flat bottom wells strips on 12 x 8 frame	10or25	100 or 200
MG04F-COOH	1x8 flat bottom wells strips on single well holding frame	10or25	100 or 200
MC04F-COOH	96 wells plate – flat bottom	10or25	100 or 200
TS 127504-COOH	12 x 75 mm tubes	100	1000-2000

BIOTIN COATED

code	description	quantity per pack	quantity per case
MT15F2-BIO	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG15F-BIO	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC15F-BIO	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127515-BIO	12 x 75 mm tubes	20	500 or 1000

CALMODULIN COATED

code	description	quantity per pack	quantity per case
MT14F2-CAL	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG14F-CAL	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC14F-CAL	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127514-CAL	12 x 75 mm tubes	20	500 or 1000

biomat

CONCANAVALIN A COATED

code	description	quantity per pack	quantity per case
MT07F2-CON A	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG07F-CON A	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC07F-CON A	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127507-CON A	12 x 75 mm tubes	20	500 or 1000

JACALIN COATED

code	description	quantity per pack	quantity per case
MT10F2-JAC	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG10F-JAC	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC10F-JAC	96 wells plate – flat bottom	1(sealed barrier bag)	100 or 200
TS 127510-JAC	12 x 75 mm tubes	20	500 or 1000

POLY-L-ARGININE COATED

code	description	quantity per pack	quantity per case
MT13F2-AR	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG13F-AR	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC13F-AR	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127513-AR	12 x 75 mm tubes	20	500 or 1000

POLY-D-LYSINE COATED (MW 70.000-150.000)

code	description	quantity per pack	quantity per case
MT17F2-LY-D	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG17F-LY-D	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC17F-LY-D	96 wells plate – flat bottom	1(sealed barrier bag)	100 or 200
TS 127517-LY-D	12 x 75 mm tubes	20	500 or 1000

POLY-L-LYSINE COATED (MW 70.000-150.000)

code	description	quantity per pack	quantity per case
MT12F2-LY-L	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG12F-LY-L	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC12F-LY-L	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127512-LY-L	12 x 75 mm tubes	20	500 or 1000

biomat

PROTEIN A COATED

code	description	quantity per pack	quantity per case
MT06F2-PA	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG06F-PA	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC06F-PA	96 wells plate – flat bottom	1(sealed barrier bag)	100 or 200
TS 127506-PA	12 x 75 mm tubes	20	500 or 1000

PROTEIN G COATED

code	description	quantity per pack	quantity per case
MT09F2-PG	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG09F-PG	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC09F-PG	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127509-PG	12 x 75 mm tubes	20	500 or 1000

STREPTAVIDIN COATED

code	description	quantity per pack	quantity per case
MT0STF2-SA5/200	1x8 flat bottom wells strips on 12 x 8 frame-coating level 200 ul*	1(sealed barrier bag)	50 or 100
MG0STF-SA5/200	1x8 flat bottom wells strips on single well holding frame- coating level 200 ul*	1(sealed barrier bag)	50 or 100
MC0STF-SA5/200	96 wells plate - flat bottom- coating level 200 ul*	1(sealed barrier bag)	50 or 100
TS 12750STF-SA5	12 x 75 mm tubes -Streptavidin coated-coating level 500 ul*	20	500 or 1000

* different coating levels are available upon request

WHEAT GERM COATED

code	description	quantity per pack	quantity per case
MT11F2-WG	1x8 flat bottom wells strips on 12 x 8 frame	1(sealed barrier bag)	50 or 100
MG11F-WG	1x8 flat bottom wells strips on single well holding frame	1(sealed barrier bag)	50 or 100
MC11F-WG	96 wells plate - flat bottom	1(sealed barrier bag)	100 or 200
TS 127511-WG	12 x 75 mm tubes	20	500 or 1000

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BIOMAT S.n.c.
 VIA ZENI 8
 I-38068 ROVERETO (TN)

È CONFORME AI REQUISITI DELLA NORMA
 HAS BEEN FOUND TO COMPLY WITH THE REQUIREMENTS OF
UNI EN ISO 9001:2000

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Alessio Galazzo
 Technical Responsible

SINCERT
INTERNATIONAL CERTIFICATION

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"The validity of the present certificate depends on the annual surveillance every 12 months and on the complete review of company's management system after three-years."

TUV Italia • Gruppo TÜV SÜD • Via Carducci 125, Pal. 23 • 20098 Sesto San Giovanni (MI) • Italia • www.tuv.it TÜV®

Biomat s.n.c. via zeni 8 I-38068 rovereto ITALY
 phone +39-0464-443320/ fax 443159
 e-mail biomat@tin.it <http://www.immunosurface.com>

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Biomat certifies that this lot has been manufactured according to applicable standard operating procedures, current ISO 9001:2000 regulation.
 This product has been tested and has passed all established specifications

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Method of analysis

Coating solution: 2.5 µg/ml of Rabbit IgG in 0.1 M Carbonate Buffer pH 9.6.
 Conjugate: Goat Anti-IgG-HRP

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Final control

This product has been tested and approved for the following physical characteristics:

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Dr. Maurizia Pettenati
 Quality Control manager

Copia conforme al certificato originale firmato e conservato presso i nostri archivi
 Original signed certificate is filed in our offices

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